

KANGDA®

Medical Brace and Support Series



Beijing Jinwei Kangda Medical Instrument Ltd.

Address: Room 1406-1409, 1st Building, No. 14 Courtyard, South
Road of Western 3rd Ring-Road, Fengtai District, Beijing 100073
Tel: +86-10-63334994 Fax: +86-10-63334994
Http://www.kdbd.com.cn Http://www.jwkd.en.alibaba.com

 Beijing Jinwei Kangda Medical Instrument Ltd.

Company Profile

We are professional manufacturer of medical polymer orthopedic instruments and sports protective equipments which specialize in these products' research and development, production, sales at home and abroad. We are the supplier of the central materials reserve, China MEHOCO CO.,LTD's export and military material procurement. We have a unity and hard work, innovation professional team, superior production and testing equipment. The company has 24 production lines of casting tape and splint, as well as 20 production lines of body support, orthosis brace. Our company cooperates with U.S Amerix company and then build R & D Center.

We have won a national patent for invention, numbers of national patents and quality excellent award of Beijing. The company has approved by ISO13485 and is in the lead of industry. All the products have CE and FDA certificates. KANGDA product has sold well for 17 years. The products are exported to USA, Europe, the Middle East, Southeast Asia and more than 40 countries and regions.

Welcome fellow members and medical organizations to co-develop and reach win-win cooperation.



Contents

Body Support Series

Knitted Support.....	02
Sport Support	04
Enhanced Support.....	06

Medical External Fixation Brace

Head & Neck Series.....	09
Thoracic Lumbar Series.....	11
Shoulder Elbow Brace Series.....	13
Wrist Brace Series.....	15
Hand Brace Series.....	16
Hip Brace Series.....	19
Knee Brace Series.....	20
Ankle Foot Brace Series.....	22



Body Support Series



Knitted Support

- Overview :**
- Select natural quality bamboo fiber, high absorption capacity, without bad smell, sweat-absorbing and cold-proof, breathable.
 - Knitted according with ergonomics, four-way-elastic, fit and comfortable.
 - Special knitted technical design matching different joints, plays the role of immobilization, protection and adjuvant therapy to joints and muscle.
 - 2 colors mix, fashion and beautiful.
 - Imported equipments, leading technology, guaranteed quality.
- Indications:**
- Applicable to the tightening and fixing damaged parts and restricts the movement of damaged parts to play the role of adjuvant therapy for partial damage.



Knee Support *KDHX-03*

Characters:

- Designed according with ergonomics, balanced pressure because of adjuvant therapy effect for knee joint diseases.
- Protect and reinforce the strength of knee joint.

Size S M L XL XXL



Waist Support *KDHY-03*

Characters:

- Improve the strength of waist with 4 internal elastic support bars, comfortable to wear and difficult to roll up.
- Adjuvant therapy effect for lumbar disorders.
- Relieve the back pain and other discomfort caused by long time work or posture.

Size S M L XL XXL



Leg Support *KDHT-03*

Characters:

- Designed according with ergonomics, balanced pressure because of "※" shape design, adjuvant therapy effect for leg diseases.
- Protect and reinforce the strength of knee joint.
- Excellent performance of heat accumulation and retention.

Size S M L XL



Ankle Support *KDHH-03*

Characters:

- Knitted according with ergonomics, balanced pressure, adjuvant therapy effect for ankle joint diseases.
- Single side edge covering technique without sewing, domestic unique, comfortable to wear and no restriction to feet.
- Protect and improve the strength of ankle joint.

Size S M L XL

Knitted Elbow Support *KDHZ-03*



Size S M L XL

Knitted Wrist Support *KDHW-03*



Size S M L

Knitted Palm Support *KDHZ-33*



Size S M L

Sport Support

- Overview :**
- Select quality spandex, nylon and other materials, four-way-elastic, comfortable and fit.
 - Knitted according with ergonomics, balanced pressure, plays the role of adjuvant therapy.
 - High temperature compression technology without sewing, comfortable to wear and no friction to skin.
 - Silica gel fits the joints tightly, plays the role of massage, fixation, shock absorption, increasing the stability joints and protection.
 - Mix of white and grey with fluorescent green points, fashion and beautiful with innervations.
 - Imported equipments, leading technology, guaranteed quality.

- Indications:**
- Applicable to the tightening and fixing damaged parts and restricts the movement of damaged parts to play the role of adjuvant therapy for partial damage.



Knee Support *KDHX-02*

Characters:

- Silicon pad fits the patella tightly and difficult to slip off, plays the role of fixation, massage, shock absorption and protection for patella.
- Metal spring support bars on both side, soft and elastic, improve the explosive power of knee joint when doing sport.
- Improve the strength of knee joint, relieve fatigue, increase endurance, protect the knee perfectly.

Size S M L XL XXL

Palm Support *KDHZ-22*

Characters:

- Disperse the pressure of joint to surrounding soft tissue, plays the role of pressure adjustment and massage, adjuvant therapy for palm.
- Release fatigue, increase endurance, protect the palm perfectly.

Size S M L





Waist Support I *KDHY-02*

Characters:

- 2 internal memory alloy support bars, plays the role of fixation, support and restriction of movement for waist.
- Show the waist circumference with front scale.
- Back cross belt design, put on and take off freely, fit the waist perfectly.
- More suitable for losing weight, exercise and protection.
- Long time of wearing is permitted.

Size S M L XL XXL



Waist Support II *KDHY-12*

Characters:

- 4 internal memory alloy support bars, plays the role of fixation, support and restriction of movement for waist more effectively.
- Suitable for lumbar disc herniation and soft tissue injury to wear.
- Long time of wearing is unpermitted.

Size S M L XL XXL



Sport Elbow Support *KDHZ-02*



Size S M L XL

Sport Wrist Support *KDHW-02*



Size S M L

Sport Ankle Support *KDHH-02*



Size S M L XL

Back Posture Orthosis

Characters:

- Made of breathable spacer and loop ect.
- Long aluminum bar support the back well.
- Easy to adjust by pulling the 8 cross pull straps in the back.

Scope of Application:

- Correct posture of people with mild humpback.
- Correct poor posture of teenagers.

Size S M L XL



Enhanced Waist Support

Enhanced Waist Support *KDHY-01*

Characters:

- Select quality fish ribbon fabric with good breathability and absorption of perspiration.
- Strong effect of fixation and protection according with ergonomics.
- 4 internal alloy support bars, restricts the movement of waist effectively and realize double strength of compression and fixation together with 2 reinforcing bands.
- Release the pressure of chest because of high waist design.
- Quality Velcro, do not scraping clothes, long service life.
- Fashion and beautiful with innervations with fluorescent green points.
- Imported equipments, leading technology, guaranteed quality.

Size S M L XL XXL




Enhanced Waist Support *KDHY-11*

Characters:

- Select quality fish ribbon fabric with good breathability.
- 5 internal elastic support bars with good effect of fixation and protection
- According with ergonomic, comfortable to wear.
- Plays the role of adjuvant therapy with 21 internal strong magnets stone.
- Humanized design with pockets to put medicine bag or heating bag.
- Imported equipments, leading technology, guaranteed quality.

Size S M L XL XXL




Magnets Stone Inside




Enhanced Knee Support *KDHX-01*



Patella Silicone Ring



Silicone Granules




Comfort Zone

Characters:


- Adopt different pressure weaving methods to increase body feeling and effectively stabilize the knee joint.
- Specially designed patella silicone ring minimizes patella decompression.
- Silica particle design, non-slip and massage effect.
- Comfort zone with special weave to eliminate stress and improve comfort.
- Highly elastic support strip provides stable support for the knee joint.

Size S M L XL XXL

Knitted Knee Support II *KDHX-13*



Special Close-edge Design



Knee Bend Design



Characters:

- Adopt different pressure weaving methods to increase body feeling.
- Special lock design, comfortable to wear.
- Adopt knee braiding technology to improve wearing comfort.

Size S M L XL XXL

Medical External Fixation Brace





Cervical Thoracic Orthosis *KDZJ-JZ-002*


Features:

- Limit cervical and high part thoracic vertebra activities.
- Head, neck and chest fixed together can effectively control cervical rotation.
- Provide stable support.
- The height adjustment of the front and back piece can be adjusted according to the patient's body shape.
- Curved support bar can provide stable support.
- Tool-free adjustment of design, easy to operate.

Scope of application:

- Cervical and high part thoracic vertebra injuries.
- Conservative treatment of cervical dislocation.
- Cervical and high thoracic vertebra injuries are fixed postoperative.

Size	Neck Circumference	Circumference From Jaw to Occiput /cm
Universal	38~40 cm	45~50



Cervical Collar *KDZJ-JZ-001*

Features:

- Limit cervical activity.
- Provide stable support.
- Throat opening for tracheotomy and care.
- Front and rear plastic support plates for increased fixing strength.
- Designed to fit the curvature of the neck, comfortable to wear.

Scope of Application:

- Neck sprain and stability fractures.
- Postoperative fixation of cervical spine.
- Chronic cervical injury.

Size	Neck Circumference/cm
S	36~38
M	38~40
L	40~42



Cervical Collar

Features:

- Well ventilated support plate to increase cool and comfort.
- Two pieces of self-adhesive support plate can adjust height freely to ensure perfect fit between chin and neck collar.
- Double buckle structure to ensure the best fixed effect.

The Scope of Application:

- Ideal use for prevention and auxiliary treatment of cervical vertebra syndrome, such as cervical spondylosis and cervical muscle strain ect and neck support for cervical spondylosis rehabilitation after surgery.

Size	Neck Circumference/cm
S	37~42
M	42~47
L	47~52
XL	52~57



Emergency Neck Collar for Pediatric *KDZJ-JZ-008*

Features:

- Limit cervical activity.
- Provide stable support.
- Maintain normal cervical alignment.
- Extra large throat opening for ease of tracheotomy and care.
- Neck support height adjustable for using in tight emergency situations.
- Simple structure and easy to operate.

Scope of Application:

- Acute cervical injury.
- Postoperative fixation of cervical spine.

Size	Neck Circumference/cm
Universal	30~36



First Aid Neck Collar (Adult) *KDZJ-JZ-009*

Features:

- Limit cervical activity.
- Provide stable support.
- Maintain normal cervical alignment.
- Extra large throat opening for ease of tracheotomy and care.
- Neck support height adjustable for using in tight emergency situations.
- Simple structure and easy to operate.

Scope of Application:

- Acute cervical injury.
- Postoperative fixation of cervical spine.

Size	Neck Circumference/cm
Universal	36~42



Adjustable Thoracolumbar Brace

Features:

- Ergonomically designed posterior plastic frame with spine relief cutout extending from inferior. Spine of scapula to sacrococcygeal junction.
- Removable plush covered anatomical design anterior panel.

Scope of Application:

- Early onset through variable stages of osteoporosis. Compression fractures of the thoracic spine.
- Kyphosis. Thoracic mechanical back pain. Fixation of pre and post surgical.

Size	Applicable Range Size/cm
XS	55-65
S	66-76
M	77-87
L	88-98
XL	99-109
XXL	110-120

Lumbar Brace KDZJ-JZ-007


Features:

- Limit lumbar flexion ,extension and rotation.
- Polymer plastic sheet, high strength, light weight.
- Horizontal size adjustable, suitable for a wide range of people.
- Two-piece front and back construction to provide a solid, stable support for the spine.
- Paste lined design, effectively increase abdominal pressure to ensure lumbar stability.
- Waist airbag design, making the patient wear more comfortable.

Scope of Application:

- Lumbar stability fracture.
- Lumbar spondylolisthesis.
- Lumbar injury fixed after surgery.
- Lumbar vertebral body lesions.

Size	Waist Circumference /cm
Universal	70-100



Medical Elastic Lumbar Corset KDZJ-JZ-006

Features:

- Limit lumbar flexion, extension and rotation.
- Built-in six removable steel bars to provide stable support for the lower back.
- Design of compression belt on the waist and back, increase abdominal pressure, reduce lumbar load.
- The main body adopts fish silk fabric, which has good air permeability.
- Ergonomic design, comfortable to wear to wear.

Scope of Application:

- Sprain on the waist.
- Chronic lumbar strain.
- Degenerative lumbar disease.
- Degenerative lumbar instability.

Size	Waist Circumference /cm
S	75-85
M	85-95
L	95-105
XL	105-115



Thoracolumbar Brace KDZJ-JZ-012

Features:

- Limit thoracolumbar flexion, extension and rotation.
- Provide stable support.
- Polymer plastic sheet, high strength, light weight.
- Horizontal and vertical size is adjustable, suitable for a wide range of people.
- Two-piece front and back construction to provide a solid, stable support for the spine.
- Paste lined design, effectively increase abdominal pressure to ensure the stability of the thoracolumbar.
- Waist airbag design, making patient wear more comfortable.

Scope of Application:

- Thoracolumbar stability fracture.
- Thoracolumbar spondylolisthesis.
- Thoracolumbar injury fixed after surgery.
- Thoracolumbar vertebral body lesions.

Size	Waist Circumference /cm
Universal	70-100





Arm Sling (Enhanced) KDZJ-ZB-004

Features:

- Secure the elbow to function.
- Reduce affected shoulder load.
- Waist band design holds the waist and forearm together.
- Prevent upper limb swing.
- Two forearm strap to increase brace stability.

Scope of Application:

- Upper extremity mild fracture fixation.
- Reset of shoulder subluxation after fixation.
- Upper limb muscle strain or sprain.
- Suspension after the plaster is fixed.

Size	Forearm Length /cm
Universal	27-40



Arm Sling (Summer) KDZJ-ZB-002

Features:

- Secure the elbow to function.
- Reduce affected shoulder load.
- Mesh material, excellent breathability.
- Simple structure, easy to wear.

Scope of Application:

- Upper extremity mild fracture fixation.
- Reset of shoulder subluxation after fixation.
- Upper limb muscle strain or sprain.
- Suspension after the plaster is fixed.

Size	Forearm Length /cm
Universal	27-40



Arm Sling (Winter) KDZJ-ZB-003

Features:

- Secure the elbow to function.
- Reduce affected shoulder load.
- Thickening warm material, suitable for winter and severe cold.
- Simple structure, easy to wear.

Scope of Application:

- Upper extremity mild fracture fixation.
- Reset of shoulder subluxation after fixation.
- Upper limb muscle strain or sprain.
- Suspension after the plaster is fixed.

Size	Forearm Length /cm
Universal	27-40

Orthopedic Arm Sling

Features:

- Sticky cloth+ sponge + sticky buckle.
- Suitable for forearm fracture and strain of wrist joint and fingers.
- It is economic, adjustable and easy to use. The patient can wear it without assistance from other people.
- Attention to keep clean in use, put in the dry and ventilated place, out of moisture.

Scope of Application:

- Shoulder dislocation.
- Shoulder joint injury postoperative fixed abduction.
- Humerus surgical neck fracture.
- Shoulder impact syndrome.
- Shoulder cuff muscle damage.



Size	Universal
Right and left are divided	

Adjustable Elbow Brace KDZJ-ZB-001

Features:

- Fix the elbow at an angle or a range of angles.
- Thickened alloy braces to provide strong support for elbow.
- Adjustable length is suitable for different arm length of patients.
- Removable hand rest can add comfort of wearing.
- Simple angle chuck, the angle adjustment is simple and convenient.

Scope of Application:

- Elbow dislocation.
- Postoperative recovery after elbow joint injury.
- Distal humerus and proximal ulnar radius were fixed.

Size	Forearm Length /cm
Universal	20-35



Angle Chuck

Adjustable Orthopedic Shoulder Support Brace


Features:

- Ok cloth + plastic ring + high quality vercle.
- Suitable for periartthritis of shoulder, and the shoulder sprain.
- Buckle closure strap for ease of application and removal.
- Attention to keep clean in use, put in the dry and ventilated place, out of moisture.

Scope of Application:

- Shoulder dislocation.
- Shoulder joint injury postoperative fixed abduction.
- Humerus surgical neck fracture.
- Shoulder impact syndrome.
- Shoulder cuff muscle damage.

Size	Universal
Right and left are divided	





Palm Wrist Brace

Features:

- Convenient and comfortable to wear.
- Pulley system strength wrist fixation.

Scope of Application:

- Trauma or postoperative of wrist chronic soft tissue injury.
- Radial nerve palsy.
- Fixation of mild or moderate wrist sprain.
- Conservative treatment of tendonitis or tenosynovitis.
- Wrist rheumatoid arthritis.
- Fixation after removing plaster bandage.

Size	Applicable Range Size /cm
S/L	Adult



Palm Wrist Orthosis KDZJ-ZW-001

Features:

- Made of high-elastic polymeric material.
- Alloy bars enclosed suitable for wrist.
- Don't restrict finger movement.
- Self-adhesive velcro strap design, easy to take off and put on.

Scope of Application:

- Mitigation and prevention for waist injury, wrist pain by over-fatigue, carpal tunnel syndrome.
- Post-operative immobilization for wrist movement control in the phase of rehabilitation.

Size	Wrist Circumference /cm
S	18-21
L	21-24



Fixed Infusion Splint KDZJ-ZW-002

Features:

- There are hardware strips designed according to the palm shape.
- Does not affect the distal phalanx free movement, improve comfort.
- Simple structure and easy to operate.

Scope of Application:

- Patients who need to fix palm wrist part when infusing.

Size	Wrist Circumference /cm
Universal	12-18



Distal Metacarpal Brace KDZJ-SB-003

Features:

- Far distal phalanx fixation.
- According to the distal phalanx shape design, comfortable to wear.
- Can be shaped to a degree based on the hand shape.
- Velcro design, fixed firmly.

Scope of Application:

- Sprain on the distal joint of finger.
- Far distal phalanx stability fracture.
- Far distal phalanx injury postoperative fixation.

Size	Phalanx Circumference /cm
Universal	4-8



Middle Metacarpal Brace KDZJ-SB-004


Features:

- Middle phalanx fixed.
- Designed according to the phalanges of the middle section, comfortable to wear.
- Can be shaped to a degree based on the hand shape.
- Velcro design, fixed firmly.

Scope of Application:

- Sprain on the middle finger.
- Middle phalanx stability fracture.
- Middle phalanx injury postoperative fixation.

Size	Phalanx Circumference /cm
Universal	4-8



Proximal Metacarpal Brace KDZJ-SB-005


Features:

- Recent phalanx fixation.
- According to the recent phalanx shape design, comfortable to wear.
- Can be shaped to a degree based on the hand shape.
- Velcro design, fixed firmly.

Scope of Application:

- Sprain of finger near section.
- Recent phalanx stability fractures.
- Recent phalanx phalanx injury postoperative fixation.

Size	Phalanx Circumference /cm
Universal	4-8





First Metacarpal Brace *KDZJ-SB-001*


Features:

- The thumb and first metacarpal bone are fixed.
- Designed for the first metacarpal shape, comfortable to wear.
- Can be shaped to a degree based on the hand shape.
- Velcro design, fixed firmly.

Scope of Application:

- Sprain of thumb and first metacarpal.
- The thumb and first metacarpal stability fracture.
- The thumb and first metacarpal are fixed postoperatively.

Size	Wrist Circumference/cm
Universal	14-24



Fifth Metacarpal Brace *KDZJ-SB-002*

Features:

- Fifth metacarpal fixation.
- Designed for the fifth metacarpal shape, comfortable to wear.
- Can be shaped to a degree based on the hand shape.
- Velcro design, fixed firmly.

Scope of Application:

- The fifth metacarpal sprain.
- Fifth metacarpal stability fracture.
- The fifth metacarpal injury is fixed postoperatively.

Size	Wrist Circumference/cm
Universal	14-24



Thumb Wrist Splint

Features:

Made of special elastic material, wrist and thumb girth is adjustable, support bar removable and can curve to fix well.

Scope of Application:

Carpometacarpal joints and metacarpophalangeal joint fracture, tendon, tendon sheath and ligament injury, inflammation caused by Strain pain.

Size	Left or Right
Free	



Finger Splint *KDZJ-SB-008*


Features:

- Finger immobilization.
- Built-in plastic support bar, moderate intensity, good support.
- The top velcro design allows for effective attachment of the strap to the hand cleat.
- Increase the fixed effect.
- Simple structure, easy to wear.

Scope of Application:

- Finger sprain.
- Fracture of finger stability.
- Finger injury is fixed after surgery.
- Adjunctive treatment of tenosynovitis.

Size	Phalanx Circumference /cm
Universal	4-8



Wearing Show of Different Parts

Trigger Finger Splint *KDZJ-SB-009*



Features:

- Finger immobilization.
- Built-in metal bracelet designed according to the shape of the metal, high strength, good support.
- Winding elastic straps for freedom of adjustment.
- Simple structure, easy to wear.

Scope of Application:

- Finger sprain.
- Fracture of finger stability.
- Finger injury is fixed after surgery.
- Adjunctive treatment of tenosynovitis.

Size	wrist circumference/cm
Universal	12-23



Wearing Show of Different Parts



Adult Hinged Hip Abduction Orthosis **KDZJ-KB-001**



Features:

- Control hip abduction, adduction, flexion, extension angle.
- Length and height is adjustable, suitable for a wide range of people.
- Swap left and right design.
- Simple structure, easy to wear.

Scope of Application:

- Conservative hip dislocation treatment.
- Hip injuries postoperatively fixed.
- Fixed after hip replacement.
- Femoral head necrosis is required for external fixators.

Size	Waist Circumference /cm
Universal	70-100



Adjustable Knee Brace **KDZJ-XB-001**

Features:

- Fix the knee joint at an angle or a range of angles.
- Tool-free adjustment of the brace length, easy to operate.
- Alloy sticks, high strength and light weight .
- Snap-in design for easy removal.

Scope of Application:

- Knee ligament injury.
- Knee stability fractures.
- Knee injury is fixed postoperatively.
- Patella fracture.

Size	Up to 15cm of Patella Center of Thigh Circumference /cm
Universal	31-70



Adjust the Chuck Precisely



Adjust Length without Tool



Two-way Interactive Magic Stick

Daily OA Knee Brace

Features:

- Provide offloading of the medial or lateral compartment for mild to moderate osteoarthritis.
- The single-sided low profile hinge encourages patients compliance.
- “U”Buttress controls the patella in a more natural position..
- Open anterior thigh section for easy donning and doffing.
- Comfortable, breathable material provides focused compression around the knee joint to address swelling and promote blood flow.

Scope of Application:

- For uni-compartment OA patient.
- Mild to moderate OA knee user.
- Help to support the knee joint and share the loading of affected joint.

Size	Thigh Circumference/cm	Calf Circumference/cm
S	30-37	25-30
M	37-44	30-35
L	44-52	35-42
XL	52-60	42-50



Children Hinged Hip Abduction Orthosis **KDZJ-KB-002**

Features:

- Control hip abduction, adduction, flexion, extension, internal rotation, external rotation angle.
- Adjustable thigh support for patients with different leg length.
- Optimized metal hinges, high strength and light weight.
- Removable liner for easy cleaning.

Scope of Application:

- Conservative treatment of pediatric developmental hip dysplasia.
- Pediatric hip dislocation after fixation.

Size	Waist Circumference /cm
S	35-45
M	45-65
L	65-70





Knee-Ankle-Foot Orthosis KDZJ-XB-006

- Features:**
- The knee flexion and extension can be adjusted according to the needs of the physician.
 - High-strength alloy struts to provide stable support.
 - Brace length, thigh length can be adjusted according to different patient leg type.
- Scope of Application:**
- External fixation of knee and tibiofibular injuries.
 - External fixation of knee and tibiofibular injuries after surgery.
 - Prevention and correction of knee contractures.

Adjustable Circumference	Adjustable Hinges
Size	Up to 15cm of Patella Center of Thigh Circumference /cm
Universal	36-57



Ligament Knee Support

- Features:**
- 1/8" thick neoprene with nylon lining.
 - Crossing straps control patella movement.
 - Open patellar design with soft-cushioned felt.
 - Upper and lower coronal fasteners offer further fit.
 - Pull-on version.
- Scope of Application:**
- Mild instabilities of knee joint.
 - Articular effusion & swelling.
 - Chronic inflammation after operation or trauma.

Size	S M L XL
------	----------



Knee Immobilizer

- Features:**
- Three piece design for precise fit.
 - Removable medial, lateral and posterior stays.
 - Comfortable, washable, skin-friendly terry toweling cloth material.
 - Velcro straps with loop closures for simple fitting.
- Scope of Application:**
- Post-operative immobilization.
 - Pre-operative immobilization due to pain.
 - Patella dislocation.

Size	12" 14" 16" 18" 20" 22" 24"
------	-----------------------------



Knee Wrap KDHXX-66

- Features:**
- 3 big straps ensure the reliable supreme support of the entire knee joint.
 - The embedded spring bars provide premium stability and shock absorption.
 - We used the best Velcro hoop and loop system which allows doing the physical activities without worrying that the brace will slip off.
- Scope of Application:**
- Protects all kinds of sport and injury, plays the role of fixing and restraining

Size	Universal (fits left or right knee)
------	-------------------------------------



Laced Shelter Ankle Brace

- Features:**
- Made of mesh fabric for great ventilation.
 - Special designed PP for better fitting, supporting and light weight.
 - Lace-up design for proper fastening as personal demand.
 - Fits right and left foot.
- Scope of Application:**
- Mild to moderate ankle distortion.
 - Chronic ankle instability.
 - Prophylactic and therapeutic use in sport activities.
 - Prevent abnormal inversion or eversion via rigid contoured uprights.

Size	S M L XL
------	----------



Adjustable Ankle Foot Orthosis KDZJ-ZH-002

- Features:**
- Unique flexion and extension adjustment band can be further adjusted according to different treatment requirements of a fixed angle.
 - The forefoot cushion provides varying degrees of traction to the achilles tendon.
 - Foot non-slip pads, effective non-slip shock absorption.
 - Bilateral design, one product fits left and right parts.
- Scope of Application:**
- Prevention of prolonged bed foot drooping.
 - Achilles tendon fixed after surgery.
 - Plantar fasciitis.

Size	Foot Length /cm
S	<19
M	20-26
L	>27



Ankle Brace with Strap

- Features:**
- Made of durable and tough nylon fabric.
 - Extra non-elastic cross-stirrup straps provide maximal stability.
 - Bilateral plastic stays for additional support.
 - Velcro lace closure design for convenient application.
- Scope of Application:**
- Mild to moderate ankle distortions.
 - Treatment of acute ankle sprain (alternative to taping)
 - Chronic ankle instability.

Size	Ankle Circumference/cm
S	15-20
M	20-25
L	25-30
XL	30-35



Foot Drop Splint KDZJ-ZH-001

- Features:**
- High resistance to bending fatigue, enabling repeated ankle dorsiflexion.
 - The rear hollow design, comfortable to wear.
 - Integrated molding, simple structure, easy to wear.
- Scope of Application:**
- Passively limit plantar flexion.
 - Simple foot drooping.
 - Feet ptosis due to common peroneal nerve injury.
 - Improvements in gait with mild hemiplegia.

Size	Shoe Size /cm
S	35-37
M	37-39
L	39-41
XL	41-44



Ankle Brace with Airbag KDZJ-ZH-006

Features:

- Thermoplastic panels provide stable ankle support.
- The medial inflator ensure fit to ankle shape to enhance the anchoring effect.
- Adjustable straps designed to be worn more applicable.
- Overlapping design to increase the comfort of brace wear.

Scope of Application:

- External fixation of ankle sprain.
- Unstable external fixation of the ankle.
- Postoperative external fixation of ankle soft tissue injury.
- External fixation after plaster removal.

Size	Ankle Circumference /cm
Universal	19-27



Adjustable Strap



Overlapping Design



Inflatable Device

Dennis Splint KDZJ-ZH-003

Features:

- Depending on the purpose of treatment to adjust foot spacing, foot valgus, foot outreach, adduction angle.
- Aviation aluminum, light weight, high strength.
- Tool-free adjustment of design, easy to operate.

Scope of Application:

- Pediatric clubfoot correction.
- Pediatric foot, valgus correction.
- Pediatric correction of the inner and outer horoscopes.

Size	Orthopedic Shoe Length /cm
8	8
10	10
12	12
14	14
16	16



Quality Creates Brand

Integrity Achieves Future